# MT110-B50A1.5-IR-Hk + MDS1C-xx

**IR TUNABLE AO MODULATOR SYSTEM** 

## Description

This complete system has been specially designed for tuneable Ti:Sa lasers. It consists of having a modulator in association with a MDS driver in order to provide constant diffraction efficiency as well as a constant output angle whatever is the laser wavelength between 690 nm and 1064 nm.

Common application can be Stimulated Emission Depletion (STED) microscopy, multi photon imaging system and many others.


## Features

- External Control/USB/RS 232
- Tuneable capabilities fro 690nm up to 1064 nm
- Constant output angle/full efficiency over whole wavelength range


Efficiency versus wavelength without input angle

readjustment

Access to your operating manual



## Frequency tracking versus wavelength



# AA OPTO-ELECTRONIC - www.aaoptoelectronic.com

AA Sa reserves the right to make changes to the products or information herein without notice. No liability is assumed as a result of their use or application. AA Sa, 18, rue Nicolas Appert F-91898 ORSAY Tel: +33 (0)1 76 91 50 12 Fax: +33 (0)1 76 91 50 31

Output Angle versus wavelength



CA



# MT110-B50A1.5-IR-HK

### **Product Overview**

These AO modulators work with a frequency range of 110 MHz +/- 25 MHz within the wavelength range 690 nm and 1064 nm. They are equipped with a passive heatsink and offer high speed application such as amplitude modulation.



## MT110-B50A1.5-IR-Hk

690 nm to 1064 nm, AR coated		
1.5 x 2 mm <sup>2</sup>		
> 85 %		
Passive Heatsink		
+10 to +40 Non condensing		
-40 to +50 Non condensing		

Options / On request			
Aperture	Wavelength	Connector	Housing

# TECHNICAL DATA SHEET 2014

AA OPTO-ELECTRONIC - www.aaoptoelectronic.com

Rise Time (Tr) is beam diameter ( $\Phi$ ) sensitive:

$$Tr = 0.66 \frac{\Phi}{V}$$

Separation angle ( $\Delta \theta$ ) is wavelength ( $\lambda$ ) sensitive:

$$\Delta \theta = \frac{\lambda F}{V}$$

Relative Efficiency / AOIVItemporal response


Amplitude modulation bandwidth ( $F_{-3dB}$ ) is rise time (Tr) sensitive:

$$\mathsf{F}_{-3\mathsf{dB}} = \frac{0.48}{Tr}$$


RF power (P) is wavelength ( $\lambda$ ) sensitive:

$$\frac{P_1}{P_2} = \frac{\lambda_1^2}{\lambda_2^2}$$

## Relative Efficiency / AOM temporal response







## AA OPTO-ELECTRONIC - www.aaoptoelectronic.com



## **Product Overview**

The MDS1C driver is based on Direct Digital Synthesizers (DDS). They produce a tuneable, stable and accurate RF frequency signal for the Ti:Sa modulator. Its design with "on the edge" technology offers unique performance in terms of accuracy, speed and stability thanks to its internal temperature correction and high linearity design.

The built in amplifier delivers the necessary RF power to drive the acousto-optic device with reduced power consumption.

The adjustment of the driver (Frequency & power ) can be done with a remote control, USB port or through the RS 232 communication.



0	
	-
	Ν

## MDS1C-B6-34-85.135 (OEM version)


Parameter	Specifications
Number of Channels	1
Frequency Range	85 MHz to 135 MHz
Frequency Stability	+/- 2 ppm/℃
Frequency Accuracy	Nom 1 KHz
Frequency Step	Nom 1 KHz
Frequency Control	Remote Control (RC 03) + USB
Power Supply	24 VDC (<0.7 A),
Rise Time/Fall time (10-90%)	< 50 ns
Modulation Input Control /External	0-10 V or 0-5V
Modulation Input Control /Internal	Remote Control (RC 03) + USB
Extinction Ratio	Nom 70 dB
Output RF power	Max 2.5 Watts
Output Impedance	50 Ω
VSWR	< 1.5/1
Input / Output Connector	DB25/SMA
Size / Weight	207 X 99 X 26.1 mm3 / 0.6 g (OEM)
Heat Exchange	Conduction through baseplate for OEM versions
Operating Temperature	10 to 40 ℃ (max Tcase 55℃)
Storage Temperature	-40 to +50 Non condensing

Options / On request	
CONTROL	RS 232
POWER SUPPLY	110-230 VAC (laboratory version)
HOUSING	Rack 19", 1U (laboratory version)
SOFTWARE	Free download on AA Website


## AA OPTO-ELECTRONIC - www.aaoptoelectronic.com








AA OPTO-ELECTRONIC - www.aaoptoelectronic.com



## OUTLINE DRAWING – MDS1C-D6-34-85.135 (Laboratory version)



AA OPTO-ELECTRONIC - www.aaoptoelectronic.com